Uslu Kobak RZ, Akkurt B. JOTCSA. 2022; 9(4): 1335-1386.

**REVIEW ARTICLE** 



### Formation and Uses of Imidazo[1,2-*a*]pyrimidines and Related Compounds: A Review Comprising Years 2000-2021

#### Rabia Zeynep Uslu Kobak\* 向, Barbaros Akkurt 🝺

Istanbul Technical University, Faculty of Science and Letters, Department of Chemistry, Maslak, 34469 Sariyer-Istanbul, Türkiye

**Abstract:** This work covers the selected synthetic papers of imidazo[1,2-*a*]pyrimidine and its derivatives between the years 2000 and 2021. Synthesis of the heterocyclic moiety, application of this scaffold to biological activities, and secondary applications like corrosion inhibition are provided. The authors hope that readers will find the treatise useful.

**Keywords:** Imidazo[1,2-*a*]pyrimidine, syntheses, biological applications, uses.

Submitted: April 29, 2022. Accepted: October 10, 2022.

**Cite this:** Uslu Kobak RZ, Akkurt B. Formation and Uses of Imidazo[1,2-a]pyrimidines and Related Compounds: A Review Comprising Years 2000-2021. JOTCSA. 2022;9(4):1335–86.

**DOI:** <u>https://doi.org/10.18596/jotcsa.1110922</u>.

\*Corresponding author. E-mail: <u>uslur@itu.edu.tr</u>.

#### INTRODUCTION

This review article has been designed in a way that the grand category, sub-categories, and examples are provided, starting from the least recent article to the most recent one.

### SYNTHETIC WORKS ASSOCIATED WITH BIOLOGICAL ACTIVITY

#### **Known Molecules**

There are known drug candidates in the market with imidazopyrimidine skeleton. Here are some examples for them (See Figure 1):

#### Uslu Kobak RZ, Akkurt B. JOTCSA. 2022; 9(4): 1335-1386.

#### **REVIEW ARTICLE**



Divaplon

Taniplon

0

Fasiplon

O







p38 kinase inhibitor

Aurora-A kinase inhibitor

Androgen receptor antagonist





for anxiety disorders

used in CK2 inhibition

Figure 1: Some known examples of imidazopyrimidines in use.

#### GABA<sub>A</sub> Ligands that are Functionally Selective

Jensen and coworkers reported that  $\gamma$ -Aminobutyric acid (GABA) is the major inhibiting neurotransmitter in the central nervous system. Imidazo[1,2-

a]pyrimidines are  $GABA_A$  receptor benzodiazepine binding site ligands. They can show functional selectivity for the a3 subtype over the a1 one (1) (See Figure 2).



**Figure 2:** Preparation of pyridylphenylimidazo[1,2-*a*]pyrimidines as GABA<sub>A</sub> receptor benzodiazepine binding ligands.

## Synthesis and Biological Activities of Some Imidazopyrimidines

An Armenian research group has studied the synthesis and biological activities of 4-substituted 6,7-dihydro-7,7-dimethyl-5-oxo-9H-pyrano-

[4',3':4,5]- thieno [3,2-e]imidazo [1,2-a] pyrimidines and their hydrochlorides and their anticonvulsant and tranquilizer properties (2) (See Figure 3).



**Figure 3:** Synthesis of pyranothienoimidazo[1,2-*a*]pyrimidine compound.

#### Antiinflammatory Response of Imidazopyrimidine Derivatives

Vidal and co-authors (2001) prepared six imidazo[1,2-a]pyrimidine (IP) derivatives and tested on leukocyte functions in vitro and inflammatory

response was also tested, this response is induced by zymosan in the mouse air pouch. The authors reported that imidazo[1,2-a]pyrimidines possess antiinflammatory potential (3) (See Figure 4).



Figure 4: Imidazopyrimidine derivatives for researching the antiinflammatory activity.

#### Synthesis of GABA Agonists

A group from Merck, Inc with Blackaby as the first author has reported the preparation of GABA-

selective agonists through a series of synthetic steps (4) (See Figure 5).



Figure 5: Preparation of tributyl Sn-C bond and reaction with aryl halide; eight instances were reported.

**Preparation of GABA**<sub>a</sub> **Agonists Against Anxiety** Jennings et al. showed that imidazo[1,2-a]pyrimidines and triazines could function as GABA<sub>A</sub> agonists for the treatment of anxiety. The group has prepared seven compounds and investigated their affinity and efficacy at a1 and a3 subtype GABA<sub>A</sub> receptors, by varying the terminal pyridine ring. The route to synthesize the molecules has been shown in the figure below (5) (See Figure 6).



Figure 6: Preparation of 3-aryl-substituted imidazopyrimidine derivative.

#### Compounds Showing Anti-inflammatory Property and COX-2 selectivity

A Chinese research group (Zhou et al.) has investigated the synthesis and anti-inflammatory activity of imidazo[1,2-a]pyrimidine derivatives positioned adjacently with two aryl groups were designed and synthesized in order to improve their anti-inflammatory activities. The synthesized compounds presented anti-inflammatory activities with some COX-2 selectivity (6) (See Figure 7).



Figure 7: Synthesis of imidazopyrimidines for anti-inflammatory activity.

#### **Preparation of New Antimicrobial Agents**

New antimicrobial agents, containing imidazo[1,2a]pyrimidine skeletons, have been synthesized by Al-Tel and Al-Qawasmeh in 2010. The antimicrobial activities were tested with Gram-positive and negative bacteria. Some of the test compounds were efficient in the inhibition of Gram-positive and -negative strains. According to the results, the substituents on phenyl rings are the determining factor of biological activity. The antibacterial activity is dramatically lowered if the molar refractivity is increased (7)(See Figure 8).

#### **Dual KSP and Aurora-A Kinase Inhibitors**

Geng and coworkers reported that four series of dihydropyrazolo[3,4-*b*]pyridines and benzo[4,5]imidazo[1,2-*a*]pyrimidines were designed

and synthesized as being KSP and Aurora-A kinase inhibitors, at the same time. This was achieved by adding some Aurora-A kinase inhibitory fragments into the KSP inhibitor body and it was named as CPUYL064. Two enzyme inhibition assays and a cytotoxicity assay in vitro were conducted with 19 target compounds. The results indicated that some target compounds could inhibit both enzymes and several compounds were significantly inhibitive against HCT116 cell line. Although they were moderate KSP and Aurora-A kinase inhibitors, two compounds displayed a remarkable cytotoxic activity in the micromolar range, especially against the HCT116 and HepG2 cell lines. Cancer treatment with dual-function inhibitors seems to be possible (8) (See Figure 3).



Figure 8: Synthesis of new antimicrobial agents containing imidazopyrimidine moiety.



**Figure 9:** Preparation of dual-function inhibitors. R<sup>1</sup> is furyl. R<sup>2</sup> is 4-methoxy, 4-nitro, 4-trifluoromethyl, 3-chloro-4-fluoro, and 4-fluoro.

### Compounds Inhibiting Wnt/b-catenin signaling pathway

In the regulation of embryonic development and tumorigenesis, Wnt/b-catenin signaling possesses an important role. When it is deregulated, severe cases like cancer occur and Wnt signaling helps establish a good platform for the pharmacologyrelated targeting of cancer. A series of imidazopyrimidines were prepared by Cosimelli ant coworkers, and they discovered that some compounds were able to inhibit the signaling pathway of Wnt/b-catenin with the help of a luciferase reporter assay. Also, cell proliferation was found in selected cancer lines and APC or b-catenin gene was mutated. The most successful compounds in the series could downregulate the expression of Wnt target genes (9) (See Figures 10 and 11).



Figure 10: Preparation of mono- and di-benzyl-substituted imidazopyrimidine compounds.



Figure 11: Benzylation of the 5-hydroxy group.

SynthesisofSomeMannichBasesIncorporating Imidazopyrimidine CoreAeluriandco-workerspreparedMannichbasescontainingimidazo[1,2-a]pyrimidineskeleton.Imidazo[1,2-a]pyrimidineskeleton.Imidazo[1,2-a]pyrimidineskeletonwaspreparedwith a one-pot,three-componentreaction,in which

а

secondary

amine

2-aminopyrimidine,

piperazine, and excess formaldehyde in methanol were used. The products were tested in vitro against three human cancer cell lines. Antiproliferative activity was found for most of the compounds. Three of them showed effective inhibition, with GI 50 values similar to the standard drug used (10) (See Figure 12).

or



Figure 12: Formation of substituted imidazopyrimidine derivatives for screening cancer cell lines.

#### Synthesis of Tricyclic Dihydropyrimidines

Some tricyclic dihydropyrimidines were synthesized by Kaur et al. using a procedure which endows a one-pot, three-component Traube-Schwarz reaction, along with  $Zn(ClO_4)_2.6H_2O$  as the catalytic entity. All the purified compounds were screened for their *in vitro* anticancer activity, against three cancer cell lines, namely prostate cancer cells (PC3), lung cancer cells (NCI-H1299), and colon cancer cells (HCT116). The authors investigated *in vitro* DNA-intercalation ability with UV-Vis spectrophotometry, and found that the compound was inserted into the DNA base pairs and interacted strongly with the DNA double helix (11) (See Figure 13).



**Figure 13:** Synthesis of tricyclic dihydropyrimidines with a three-component reaction for the hopes of anticancer activity.

### Antigiardial Response of Imidazopyrimidines

Velázquez-Olvera and co-workers investigated in vitro of the antigiardial response of a group of imidazopyrimidine compounds against *Giardia lamblia* WB's trophozoites and albendazole was employed as the reference drug. The synergism of albendazole with the most active imidazopyrimidine compound was also studied. The most potent compound had 3-hydroxymethyl-4-fluorophenyl substituents on the skeleton, and the parent compound to compare had p-tolyl substituent only. As a conclusion, researchers found that 2-aryl-3-hydroxymethyl substituents created a good synergism against the reference (12) (See Figure 14).



**Figure 14:** 2-Aryl-3-hydroxymethyl-substituted imidazopyrimidine compound and albendazole as a reference compound for antigiardial response.

#### Synthesis and Biological Activity of Biphenyl-Containing Imidazopyrimidines

Al-Lami and coworkers reported the four-step synthesis of fourteen compounds. In the first step, 2-aminopyrimidine and biphenyl phenacyl bromide were reacted to give the imidazopyrimidine skeleton. In the second step, Vilsmeier-Haack reaction was employed to introduce an aldehyde group at the 3-position. The aldehydes formed were reacted with various aromatic amines to form

Schiff's bases in the third reaction. Sodium borohydride was used to reduce the Schiff's bases to 3-aminomethyl-2-biphenylimidazo[1,2-

*a*]pyrimidines. Different bacteria were tested for antimicrobial effect of some of the synthesized compounds (13) (See Figure 15).



Figure 15: Four step synthesis of fourteen compounds for antimicrobial assay.

#### **PDE10A** enzyme inhibitors

Moszczyński-Pętkowski and coworkers synthesized new substances containing imidazo[1,2-*a*]pyrimidine skeleton for possible interaction with PDE10A. Some

compounds displayed a very high activity and also presented a good metabolic stability (14) (See Figure 16).



**Figure 16:** Synthesis of an imidazopyrimidine compound having a triazolopyrimidine substituent and another heterocycle at 2 position for PDE10A inhibition.

#### Antitumor Activity along with Antimicrobial and Antioxidant Effect of some Imidazopyrimidine Derivatives

Farag and coworkers investigated the potential utility of imidazo[1,2-*a*]pyrimidine derivatives. Many compounds synthesized in this study showed a very

good *in vitro* antitumor activity against MCF-7 cell line. Their antimicrobial and antioxidant activities were also found to be high. Some computational studies like DFT and SAR were also reported (15) (See Figures 17).





Figure 17: Synthesis of benzimidazopyrimidines from nitriles and aminobenzimidazoles.

### Homopiperazine-Containing

**Imidazopyrimidine Derivatives** Homopiperazine-containing new imidazo[1,2*a*]pyrimidine derivatives (3 in total) were prepared and screened by Mantipally and coworkers, in terms of *in vitro* cytotoxicity against HeLa and A549 by MTT technique (16) (See Figure 18).



**Figure 18:** Preparation of imidazopyrimidine from aminoimidazole and substitution of the chloro group with homopiperazine and further functionalization of the homopiperazine side group.

#### Dihydropyrimidine as Monastrol Analogs

Tawfik and co-workers prepared twenty-one 3,4dihydropyrimidine compounds, as monastrol analogs, having 1,3-benzodioxole moieties at position 4 and also having different substituents at positions 2, 3 and 5. These compounds were screened towards 60 cancer cell lines under NCI (USA) protocol. The researchers also reported the assay of levels of active caspase-3 and caspase-9. They also studied molecular docking of some compounds to illustrate the interaction between inhibitors and the kinesin spindle protein allosteric binding site (17) (See Figure 19).



Figure 19: Synthesis of imidazopyrimidine derivatives as monastrol analogs.

Synthesis and Antimicrobial Activity of Quinoline-Containing Imidazopyrimidines Quaternary salts from a quinoline compound were

prapared by Elenich and coworkers, then imidazo[1,2-a]pyrimidine derivatives, along with other heterocycles, were produced with the use of dinucleophilic reagents. A great deal of the synthesized compounds showed antimicrobial activity (18) (See Figure 20).



Figure 20: Imidazopyrimidine compounds for antimicrobial activity.

#### Fluoromethylated Imidazopyrimidines

Jismy and coworkers synthesized fluoromethylated derivatives of imidazo[1,2-*a*]pyrimidines and benzimidazo[1,2-*a*]pyrimidines with a Michael addition/intramolecular cyclization. In this reaction, 2-aminoimidazole derivatives were reacted with 4,4,4-trifluorobut-2-ynoate under C-O bond

activation. The formed compounds are examples of C-C, C-N, and C-S bond-forming synthons. Some derivatives possessed micromolar IC50 values against monoamine oxidase B and butyrylcholinesterase and as known well, these two enzymes are important targets considered by neurodegenerative disorders (19) (See Figure 21).

Uslu Kobak RZ, Akkurt B. JOTCSA. 2022; 9(4): 1335-1386.

#### **REVIEW ARTICLE**



Figure 21: Fluoromethylated derivatives of imidazopyrimidines.

#### Stabilizing Effect for KRAS 24

Pancreatic, pulmonary, and colorectal cancers involve an oncogene named KRAS, so it constitutes an important target for the treatment of these cancers. A group of molecules were found by D'Aria and coworkers to perform a stabilizing effect on the NMR structure of 22-mer KRAS 24. Finally, compound 19 was found to be the best candidate and it may be considered that this compound is possibly a prototype of a new class of drugs which could be used for the treatment of tumors which express mutations of KRAS and which resist current therapies (20) (See Figure 22).



Figure 22: Quaternarization and ring formation from a benzimidazopyrimidine.

### Ten Year Review of Therapeutic Potentials of Pyrimidines

Zhuang and coworkers described the 10-year review of therapeutic potentials of pyrimidine-based

compounds and their antimicrobial activities (21) (See Figure 23).



Figure 23: Antimicrobial imidazopyrimidines reviewed by Zhuang and coworkers.

#### **HIV-1** Integrase Inhibitors

As known well, anti-human immunodeficiency virus 1 (anti-HIV-1) drugs make it a target for HIV-1 integrase (IN) enzyme. The researchers with Zadeh as the first author modified the structure of N-

arylindole  $\beta$ -diketoacids (favorite inhibitors of IN) with 4-oxo-4,10-dihydrobenzo[4,5]imidazo[1,2a]pyrimidine-3-carboxylic acid derivatives for the development of new anti-HIV-1 drugs (22) (See Figure 24).



Figure 24: Synthesis of benzimidazopyrimidine compounds for HIV-1 integrase inhibitors.

### Attempts to Use Imidazo[1,2-*a*]pyrimidines in Photodynamic Therapy

Some imidazo[1,2-a]pyrimidines were found by Lima and coworkers to generate singlet oxygen and thus, they could serve as intracellular photosensitizers. Photodynamic therapy of these compounds indicated that they absorbed and emitted within 400-500 nm, they showed low cytotoxicity when not in use (in the dark), they are taken inside cells efficiently, they fluoresce inside cells, they generate singlet oxygen with irradiation of the system, and cancer cells are killed in two hours and at a low concentration (23) (See Figure 25).

#### **PURELY SYNTHETIC WORKS**

#### Vilsmeier Formylation of Imidazo[1,2a]pyrimidines

Vilsmeier reaction is largely employed to add formyl groups to a molecule. In this reaction, phosphorus oxytrichloride and N,N-dimethylformamide are used. A Latvian group with Saldabol was the first author has used this methodology to add formyl (-CH=O) groups to imidazo[1,2-*a*]pyrimidine at the 3

position. When the formylating reagent is taken as double, the side furyl group was also formylated (24) (See Figure 26).

#### Synthesis of 2-carbamoyl-3-arylimidazo[1,2a]pyrimidines

Majcen-Le Marechal and coworkers have reported a series of 3-aryl-substituted imidazo[1,2a]pyrimidine-2-carboxamide structures (25) (See Figure 27).

#### Synthesis of Dihydroimidazopyrimidine

Kochergin coworkers 2,3and used dihydroimidazo[1,2-a]pyrimidine in the synthesis where 2-aminopyrimidine was used by two different routes. The first route involved reacting 2aminopyrimidine with ethylene halohydrin (bromo or chloro) and further reacting with thionyl chloride, while the other method involved 1,2-dibromoethane reacting with 2-aminopyrimidine and when the intermediary product was treated with sodium hydroxide and then with hydrobromic acid, dihydroimidazopyrimidine was obtained (26) (See Figure 28).



Figure 25: Synthesis of a tricyclic hetero compound for photodynamic therapy applications.



Figure 26: Formylation of 2-furyl-containing imidazopyrimidines.



Figure 27: A series of 3-aryl-substituted imidazopyrimidine-2-carboxamides.



Figure 28: Synthesis of 2,3-dihydroimidazo[1,2-a]pyrimidine with two different routes.

#### Fusion to the imidazopyrimidine core

A double annelation reaction by Chowdhury and colleagues yielded a variety of tri- and tetracyclic heterocycles in a one-pot reaction (27) (See Figure 29).

### Antipyrine-containing imidazopyrimidine heterocycle

Abu-Elmaati described the syntheses of many heterocycles, in which he reported an imidazo[1,2-*a*]pyrimidine that contains an antipyrine substituent (28) (See Figure 30).

#### One-pot synthesis of imidazo[1,2a]pyrimidines

Late professor Alan Roy Katritzky is a legend in heterocyclic chemistry, and he has reported the one-pot synthesis of imidazo[1,2-a]pyrimidines. He and his colleagues used 2-aminopyrimidine and 1,2bis(morpholinyl)-1,2-dibromoethane, with which two intermediates, namely A and B occurred with the sequential removal of benzotriazolyl (denoted as Bt here) moieties, then a rearrangement in the ring gives two imidazopyrimidine derivatives (29) (See Figure 31).



Figure 29: Tricyclic and tetracyclic imidazopyrimidine derivatives in one-pot reactions.



Figure 30: Synthesis of pyrazole-containing imidazopyrimidines.



Figure 31: Syntheses of two imidazopyrimidines with the use of dibromodimorpholinylethane. Two compounds were reported.



Figure 32: Synthesis of imidazopyrimidines on solid support.

Heterocycle formation on a solid support

El Kazzouli and his coworkers reported the formation of imidazo[1,2-a]pyrimidine by employing a solid-phase support (30) (See Figure 32).

#### Formylation on Monosaccharides

Bari et al. reported another formylation reaction applied on monosaccharides and the heterocycle was obtained with 2-aminobenzimidazole (31) (See Figure 33).



Figure 33: Imidazopyrimidine obtained with monosaccharides.

#### Imidazo[1,2-a]pyrimidin-5-ones

According to the publication by Sączewski et al., imidazolin-2-oxime ethers react with alkenes and imidazo[1,2-*a*]pyrimidin-5-ones. In addition, a retro-ene reaction takes place when ethyl prop-2ynoate and 2,3-dihydroimidazo[1,2-a]pyrimidin-5(1H)-one is obtained (32) (See Figure 34).

#### Fluorescent Probes for the Localization and Function of Peripheral Benzodiazepine Receptor (PBR)

PBRr is mainly found in the mitochondria, and its overexpression is encountered, for example, in glioma, breast cancer, Alzheimer's disease, and activated microglia. Among the handful of compounds synthesized by Laquintana and coworkers, an imidazopyridine-7-nitrofurazan conjugate showed the best result and could be considered as a new useful fluorescent probe for the visualization of activated microglia and PBR (33) (See Figure 35).

### Reaction of an Imidazole and a Pyran to Give Imidazopyrimidine Derivatives

A new method was proposed by El Otmani and their co-workers used 4-hydroxy-6-methylpyran-2-one and 2-amino-4,5-dicyanoimidazole with an alcohol (methanol, ethanol, *n*-propanol, and *n*-butanol) to obtain dicyanoimidazopyrimidine acetic acid derivatives and their corresponding esters of the alcohol used (34) (See Figure 36).



Figure 34: Use of imidazolin-2-oxime ethers in the construction of partly hydrogenated imidazopyrimidines.

### Formation of Imidazopyrimidines via a Nucleophilic Attack on Styrenes

A new stereoselective nucleophilic attack on styrenes was described by Muzalevskiy and colleagues. First, the styrene is reacted with sodium methoxide and potassium tert-butoxide, then brominated to yield the trifluoromethyl ketones. These compounds were found to be good starting materials having a trifluoromethyl moiety. Imidazopyrimidine compounds were reported to be obtained in medium yields (35)(See Figure 37).

#### Hydrazinolysis of Imidazopyrimidines

Ermolat'ev and co-workers reported the microwaveassisted hydrazinolysis of imidazopyrimidines to yield mono- and di-substituted imidazoles. This protocol is advantageous due to the absence of the use of strong acids and is better than the classical reaction in which a-haloketones react with N-acetylguanidine (36)(See Figure 38).

#### Synthesis of Benzylimidazopyrimidines

Bakherad devised an efficient synthesis of 2benzylimidazo[1,2-*a*]pyrimidines in water, incorporating Sonogashira coupling with various aryl iodides, and without copper co-catalyst, in the presence of  $K_2CO_3$  as the base. (37) (See Figure 39).



Figure 35: Synthesis of fluorescent probes for peripheral benzodiazepine receptor.



Figure 36: Formation of 7-methyl-2,3-dicyanoimidazo[1,2-a]pyrimidine-5-acetic acid and its ester.



Figure 37: Use of styrene derivatives to yield imidazopyrimidines.



Figure 38: Hydrazinolysis of imidazopyrimidines.



Figure 39: Sonogashira-powered cyclization of 2-aryl substituted imidazopyrimidines.

## Catalytic Chalcogenation of Imidazopyrimidines

Li and coworkers reported a very potent and environmentally friendly catalytic chalcogenation of imidazopyrimidine compounds with dichalcogenides. They used copper(I) iodide as catalyst and the reaction was conducted under air. Smooth reactions were observed and the yields were moderate to excellent and no other additive was required (38) (See Figure 40).



Figure 40: Chalcogenation at the 3-position by using diphenyldisulfide and copper iodide.

# Synthesisof4-Amino-1,2-dihydrobenzo[4,5]imidazo[1,2-a]pyrimidine-6-carbonitrile Derivatives

Hu and coworkers developed a fast and green method for the preparation of 4-amino-1,2dihydrobenzo[4,5]imidazo[1,2-*a*]pyrimidine-6carbonitrile derivatives. The reaction is characterized by its one-pot condensation of 2aminobenzimidazole, several aldehydes, and malononitrile and the catalyst employed was ammonium acetate in ethanol. The catalyst is very cheap and is readily available, and the workup is environmentally friendly (39) (See Figure 41).



**Figure 41:** Three-component reaction of aminobenzimidazole, an aldehyde, and malononitrile in ammonium acetate medium. R<sup>1</sup> is phenyl, 4-methylphenyl, 4-methoxyphenyl, 4-chlorophenyl, 3-chlorophenyl, 2-chlorophenyl, 4-hydroxyphenyl, 4-cyanophenyl, and 4-pyridyl.

### Synthesis of Some Imidazopyrimidines and Catalytic Hydrogenation of them

Some new derivatives of imidazo[1,2-a] pyrimidine have been studied by Borisov et al in 2013, and

catalytic hydrogenation methods were developed for them. For the first time, partial reduction of the imidazole ring on imidazo[1,2-*a*]pyrimidine hydrogenation was reported (40) (See Figure 42).



Figure 42: Formation of some imidazopyrimidine heterocycles and partial hydrogenation of some of them.

#### Synthesis of New Benzimidazole Derivatives

Novel benzimidazole derivatives were synthesized by Mehranpour and Zahiri with the reaction between 2-substituted 1,3-bis(dimethylamino)-trimethinium salts and 2-aminobenzimidazole in the presence of acetic acid or triethylamine in acetonitrile as solvent (41) (See Figure 43).



Figure 43: Formation of imidazopyrimidines via dimethylaminomethinium salts.

#### Thermal Cyclization of Thioxopyrimidines

Burbuliene et al. reported the reaction of methyl (2methylsulfanyl-6-phenyl-4-thioxopyrimidin-3(4H)yl)acetate with hydrazine hydrate, at room temperature, in 1-butanol as solvent to yield the hydrazide. They committed the same reaction at reflux temperature by employing different solvents like methanol, acetonitrile, n-butanol or dimethylformamide and they found that 1-amino-7-phenyl-5-thioxo-1,5-dihydroimidazo[1,2-*a*]pyrimi-din-2(3H)-one with a thermal cyclization pattern (42) (See Figure 44).



Figure 44: Synthesis of imidazopyrimidines with thermal cyclization of thioxopyrimidines.

#### **Review about C-H Arylation**

Aziz and Piguel, in their mini-review published in 2017, studied C-H arylations of heterocycles, in

which imidazo[1,2-*a*]pyrimidines played a major role. Some examples were provided (43). (See Figures 45-54).

Het(Ar)Br (1.15 equiv.)  
Pd(OAc)<sub>2</sub> (2 mol%)  
PPh<sub>3</sub> (4 mol%)  
Cs<sub>2</sub>CO<sub>3</sub> (2 equiv.)/dioxane or  

$$K_2CO_3$$
 (2 equiv.)/DMF  
HBr  
HBr  
HBr













Figure 48: Heteroannulation at 3- and 5-positions.



Figure 49: Development of a new strategy over the problematic, previous one.



Figure 50: Dimroth rearrangement and subsequent C-H arylation of an imidazopyrimidine.



Figure 51: Arylation and subsequent hydrazinolysis of an imidazopyrimidine.



Figure 52: C-H (het)arylation and subsequent hydrazinolysis.



Figure 53: Arylation and functionalization of a pyrimidinyl imidazopyrimidine.



**Figure 54:** Addition of a carbonyl group at 3-position of imidazopyrimidine in the presence of iron(III) chloride and aerial oxygen.

**Formation of Some Bridgehead Heterocycles** Reacting 2-aminopyrimidine and acetyl acetone, Bhagat and Telvekar reported the formation of Imidazo[1,2-*a*]pyrimidines at around 70% reaction yield (44) (See Figure 55).



Figure 55: Acetyl acetone and ethyl acetoacetate were used to prepare imidazopyrimidines.

## Use of Copper Oxide Nanoparticles in the Formation of Imidazopyrimidines

Rawat and Rawat studied the copper oxide nanoparticles-catalyzed reaction employed 2-aminobenzimidazole, several aldehydes, and several

terminal alkynes to give imidazo[1,2-*a*]pyrimidines. The reaction did not employ solvents, and 6-endodig cyclization occurs. Of the 20 compounds, 2,4bis(4-methoxyphenyl)-substituted compound was found to sense zinc ion in a fluorescent manner. WHO's maximum allowable zinc concentration in

drinking water is much higher than the detection limit of this sensor (45) (See Figure 56).



Figure 56: A three-component reaction catalyzed by HS-CuO to yield imidazopyrimidines.

C-H Bond Chalcogenation of Imidazopyrimidines Bettanin and coworkers studied the C-H chalcogenation of N-heteroaryl compounds with ammonium iodide as the catalyst, in a minimal amount of dimethylsulfoxide, water, and acetic acid (2.5 M, 2.5 M, 1 M, respectively) and no metallic catalyst was used. This approach was also very efficient in the C-H bond chalcogenation of other 5-membered N-heteroaryl compounds (46) (See Figure 57).



Figure 57: C-H chalcogenation of imidazopyrimidines.

#### Chromen-2-one-Containing Imidazopyrimidine Derivatives

Yagodinets et al. Reported the reactions of 4-(4bromoacetylphenyl)-3-hydroxy-2H-chromen-2-one with suitable heterocycles yielded quaternary salts, and the reactions of the same bromo derivative with pyrimidin-2-amine provided corresponding imidazo[1,2-*a*]pyrimidine derivatives (47) (See Figure 58).



**Figure 58:** Formation of imidazopyrimidine from a chromen-containing compound.

### IBX/NIS-Induced Intramolecular Oxidative Annulation without Metallic Catalyst

Imidazo-heterocycles were synthesized by Makra and coworkers, without a metallic catalyst, with an IBX/NIS-induced intramolecular oxidative annulation of Mannich substrates. The reaction includes iodination, NH-oxidation, formation of intramolecular C-N bonds, and retro-Claisen-Schmidt reaction with high reaction yields. A sequential one-pot reaction process has also been reported (48) (See Figure 59 and 60).



Figure 59: A three-component one-pot reaction involving the formation of an imidazopyrimidine.



Figure 60: IBX- and NIS-mediated reaction yielding imidazopyrimidine derivatives.

A Review About Pyrimidine-Containing Five-Membered Heterocycles

Pyrimidine-containing five membered heterocyclic ring systems were reviewed for the last two decades by Maji in 2019 (49). (See Figures 61-65).



Figure 61: Formation of imidazopyrimidines via an iminoalkene.



Figure 62: Using fluoro or difluoropyridine, the authors obtained imidazopyrimidines at room temperature.



Figure 63: Three-component syntheses in which magnetite was used as a catalyst.



Figure 64: Three-component synthesis of a benzimidazopyrimidine ketone.



**Figure 65:** Synthesis of benzimidazopyrimidine derivatives from an aldehyde and an acetylene derivative, with CuO nanoparticles as the catalyst.

#### **Improved Synthesis of 18F-GTP1**

In the past, 18F-GTP1's synthesis was hard and yields were low, chromatographic separations were inefficient, and product quality was variable. This synthesis was referred to as the first-generation. White and coworkers devised a more successful, second-generation synthesis, which could be extended to the kilogram scale. Lithium aluminum deuteride caused the placement of geminal deuterium atoms on the structure, then an efficient amide-forming reaction led to the key acrylamide coupling partner compound. Tricyclic imidazo[1,2a]pyrimidine compound was introduced with a highly successful annulation reaction. Unlike the first-generation synthesis, no chromatographic separations were required and the final step involved a successful and reproducible tosylation step (50) (See Figure 66).

#### Dichalcogenation of Imidazopyrimidines

Obah Kosso and coworkers treated imidazo[1,2a]pyrimidine and similar heterocycles with dichalcogenation at C6 position and the oxidizing medium was iodine coupled with dimethylsulfoxide and the starting compound were diarylchalcogenides. This strategy was efficient, and the regioselectivity and yields were excellent (51) (See Figure 67).

#### Halogenation of Imidazole-Derived Heterocycles

Neto and coworkers reported trihaloisocyanuric acid as a practical and eco-friendly reagent for regioselectively halogenation of imidazole-derived heterocycles (52) (See Figures 68 and 69).

#### **Review of Imidazopyrimidines for Direct C-H Bond Functionalization**

Patel and coworkers prepared a review article in which direct C-H bond functionalization of imidazopyrimidines and other related compounds (53) (See Figures 70-77).



Figure 66: Improved synthesis of an imidazopyrimidine derivative.



Figure 67: Dichalcogenation of imidazopyrimidine at 3-position and 6-position.



Figure 68: Use of trihaloisocyanuric acid with a heteroaromatic compound for imidazole-derived heterocycles.



Figure 69: Synthesized imidazopyrimidine bases.



Figure 70: C-H bond functionalization of an imidazopyrimidine derivative.



Figure 71: Substitution of 3-H with a hydrocarbon chain.



Figure 72: Functionalization of 3-H with an aminothiazole derivative.



**Figure 73:** Functionalization of 3-H position with a photochemical route employing an iridium complex.

![](_page_35_Figure_4.jpeg)

Figure 74: Electrolysis of an imidazopyrimidine derivative, yielding a 3,3'-dimer.

![](_page_35_Figure_6.jpeg)

Figure 75: Addition of an aminoheterocycle at 3 position with tert-butylhydroperoxide.

![](_page_35_Figure_8.jpeg)

Figure 76: Addition of an aldehyde group to the 3-position.

![](_page_35_Figure_10.jpeg)

Figure 77: Functionalization of imidazopyrimidine 3-position with a sulfonimidoyl group.

**Formation of Benzoimidazotriazolopyrimidines** Fedotov and coworkers described a new synthetic methodology to benzo[4,5]imidazo[1,2-*a*] [1,2,3]triazolo[4,5-*e*]pyrimidines. The synthesis included the condensation of 3-(arylazo)benzo[4,5]imidazo[1,2-*a*]pyrimidine-4amines from aminobenzimidazoles and 3-oxo-2phenylazopropionitrile and, with copper(II) acetate as catalyst, the mentioned compounds formed with oxidative cyclization with around 75% yield. The compounds showed good quantum yields, maximal absorption within 380-400 nm, and maximal emission within 470-500 nm (54) (See Figure 51).

![](_page_36_Figure_2.jpeg)

Figure 78: Synthesis of a benzoimidazotriazolopyrimidine compound.

#### ONE-POT, MULTICOMPONENT SYNTHETIC METHODS

#### **Multicomponent Reaction Approach**

A one-pot approach to substituted imidazo[1,2a]pyrimidines is described by Kiselyov and Smith, II, with a convenient protocol. The reaction is general with nitrile, aldehyde, and aminoheterocycle components. The reaction yields are around 70% and isolation is easy (55) (See Figure 79).

![](_page_36_Figure_8.jpeg)

Figure 79: Structures of the formed compounds.

#### A Multicomponent Reaction with Isonitriles

2-Aminopyrimidine, aldehydes, and isonitriles undergo a multicomponent reaction to yield imidazopyrimidines. A Russian research group, with Parchinsky as the first author, published a paper in which they showed the very possibility of the formation of several isomeric compounds (56) (See Figure 80).

#### One-pot multicomponent synthesis of 3amino-2-arylimidazo[1,2-a]pyrimidine derivatives

Adib and co-workers used 2-aminopyrimidine, a benzaldehyde (substituted or not), and imidazolin-2,4,5-trione in a one-pot, solvent-free and multicomponent reaction to obtain imidazopyrimidine derivatives (57) (See Figure 81).

![](_page_36_Figure_14.jpeg)

Figure 80: A multi-component reaction to give imidazopyrimidine derivatives.

![](_page_37_Figure_2.jpeg)

Figure 81: Formation of 2,3-disubstituted imidazopyrimidines via two different pathways.

#### Synthesis of aryl-benzo[4,5]imidazo[1,2a]pyrimidine amines

Reddy et al. grafted tetraethylene glycol-bridged 1vinylimidazolium mesylate onto the surface of crosslinked polyethylene glycol methacrylate. This structure served as an excellent catalyst for the synthesis of N-methyl-2-nitro-arylbenzo[4,5]imidazo[1,2-*a*]pyrimidine amines in a multi-component reaction employing 1H- benzo[d]imidazol-2-amine and (E)-N-methyl-1-(methylthio)-2- nitroethenamine and several aldehydes under solvent-free conditions. The catalyst can be reused and its catalytic activity can maintain its superiority after seven cycles. The whole protocol is green, no waste side products, no solvent involved, and no column chromatography is needed (58) (See Figure 82).

![](_page_37_Figure_7.jpeg)

Figure 82: Synthesis of N-methyl-2-nitroarylbenzo[4,5]imidazo[1,2-a] pyrimidine amines.

Synthesis of Benzoimidazopyrimidine Derivatives with NiFe<sub>2</sub>O<sub>4</sub> as the Catalytic Entity Hamidinasab and coworkers prepared a modified NiFe<sub>2</sub>O<sub>4</sub> nanoparticle for use in the multicomponent synthesis of benzo[4,5]imidazo[1,2-a]pyrimidine derivatives under green conditions. The reaction provided high reaction yields, reaction times were shorter, purification was simple, and reaction conditions were mild by environmental means (59) (See Figure 83).

![](_page_38_Figure_2.jpeg)

Figure 83: Synthesis of benzo[4,5]imidazo[1,2-a]pyrimidine derivatives.

#### Silica Sulfuric Acid as the Catalytic Entity

Basyouni and coworkers devised a simple and efficient, and eco-friendly synthetic method for synthesizing benzo[4,5]imidazo[1,2-*a*]pyrimidine-3-carbonitrile derivatives with excellent yields. They prepared the final compounds with a multicomponent reaction among 2-

aminobenzimidazole, several aldehydes, malononitrile or ethyl cyanoacetate as active nitriles, and silica sulfuric acid/ethylene glycol as the catalytic system. The significant aspects of this synthetic methodology are very short reaction times in terms of minutes, high yields, and costeffectiveness (60) (See Figure 84).

![](_page_38_Figure_7.jpeg)

Figure 84: Synthesis of imidazopyrimidine carbonitrile derivatives in two different routes.

### $Fe_3O_4@C-SO_3H$ Nanoparticles as the Catalytic Entity

In order to synthesize a separable magnetic  $Fe_3O_4@C-SO_3H$  nanoparticle (MNPs) with high catalytic activity in organic transformation, three environmental-benign and low-cost sulfonic acid

functionalized magnetic nanoparticles (Fe<sub>3</sub>O<sub>4</sub>@C-SO<sub>3</sub>H) were successfully synthesized by Damghani et al. Benzo[4, 5] imidazo[1,2-a]-pyrimidine derivatives were synthesized under solvent-free conditions in excellent yields (61) (See Figure 85 and 86).

![](_page_38_Figure_12.jpeg)

Figure 85: Magnetite-catalyzed synthesis of imidazopyrimidine compounds.

column

![](_page_39_Figure_2.jpeg)

Figure 86: An examplary reaction involving the three-component synthesis.

#### Use of Starch-Functionalized Magnetite nanoparticles were used as the catalytic entity in Nanoparticles as the Catalytic Entity the reaction. The developed reaction system used mild reaction conditions, isolation was easy, atom multicomponent synthesis One-pot, of imidazopyrimidine compounds was realized in an economy was high, products were obtained in good efficient and environmentally friendly fashion by excellent reaction to Verma et al. The components were aromatic chromatography was involved, and the catalyst was aldehydes, active methylene compounds, and 2magnetically separable and reusable. (62) (See aminobenzimidazole and ultrasonic irradiation was Figure 87). Starch-functionalized employed. magnetite starch-functionalized $Fe_3O_4$ (4 mg)

![](_page_39_Figure_5.jpeg)

![](_page_39_Picture_6.jpeg)

vields,

no

Figure 87: An aldehyde, an active methylene compound, and benzimidazole react to give a benzimidazolo[1.2-a]pyrimidine compound at room temperature, with silicafunctionalized magnetite particles as the catalytic entity. 21 compounds were reported.

#### Synthesis of Imidazopyrimidine-Based Pyrans

Imidazo[1,2-a]pyrimidine-2-carbaldehyde,

malononitrile, acids in which C-H activition is possible and sodium carbonate as base were reacted by Güngör in a one-pot fashion to yield new imidazo[1,2-a]pyrimidine-based pyrans at room temperature and similar to other compounds, a

multicomponent reaction scheme was adopted. Mild reaction conditions, the absence of the need for column chromatographic separations, and the medium to good reaction yields are the main advantages for this multicomponent reaction (63) (Figure 88).

![](_page_39_Figure_12.jpeg)

Figure 88: New imidazopyrimidine-based pyran compounds.

Heteroannulation of Imidazopyrimidines

Akbari and coworkers prepared a nanomaterial heterogeneous catalyst and used this in the Uslu Kobak RZ, Akkurt B. JOTCSA. 2022; 9(4): 1335-1386.

#### **REVIEW ARTICLE**

heteroannulation of imidazopyrimidines with a multicomponent reaction among 2aminobenzimidazole, compounds having acidic C-H entities, and several aromatic aldehydes. The protocol was green and catalyst loading was low, yields were high to quantitative, reaction times were short, and the catalyst could be recycled at least four runs (64). (See Figure 89).

#### Dihydrobenzoimidazopyrimidine **Synthesis** with ZnFe<sub>2</sub>O<sub>4</sub> as the Catalytic Entity

2-Amino-4-substituted 1,4dihydrobenzo[4,5]imidazo[1,2-a]pyrimidine-3carbonitriles were synthesized with а multicomponent reaction among benzimidazole, aldehydes, and malononitrile under ultrasonic irradiation and heterogeneous ZnFe<sub>2</sub>O<sub>4</sub> as the synthetic nanocatalyst. The methodology is especially remarkable in that the reaction time is short, handling is easy, simple, efficient, and the reaction yield is high and the catalyst is recoverable (65) (See Figure 90).

![](_page_40_Figure_5.jpeg)

Figure 89: Three-component synthesis of heteroannulated imidazopyrimidines.

![](_page_40_Figure_7.jpeg)

![](_page_40_Picture_8.jpeg)

![](_page_40_Picture_9.jpeg)

reactor, 70 °C

Figure 90: 2-Amino-4-substituted 1,4-dihydrobenzo[4,5]imidazo[1,2-a]pyrimidine-3carbonitrile compounds.

#### Magnetite@Titanium Dioxide@Ionic Liquid Zeolitic Imidazolate Framework as the **Catalytic Entity**

One-pot synthesis of benzoimidazopyrimidine compounds was utilized with a magnetite@titanium dioxide@ionic liquid-zeolitic imidazolate framework, which showed high catalytic activities. The threecomponent reaction (various aldehydes, 2aminobenzimidazole, and ethyl acetoacetate) was

conducted at room temperature. Similarly to other magnetic catalysts, the nanocatalyst was easily recovered with a magnet and the catalytic behavior of it remained the same for six catalytic runs. The advantages were that yields were excellent, procedures were simple, reaction times were short, workup was simple, and reaction conditions were mild (66) (See Figures 100 and 101).

![](_page_41_Figure_2.jpeg)

Figure 91: A part of the catalyst used visualized.

![](_page_41_Figure_4.jpeg)

Figure 92: The three-component reaction and the product obtained.

Synthesis of Imidazopyrimidines with a Ferrocene-containing Zinc(II) Dithiocarbamate complex as heterogeneous bifunctional catalyst Anamika and coworkers prepared a ferrocenedithiocarbamate zinc(II) complex to serve as a heterogeneous bifunctional catalyst in the threecomponent synthesis of imidazopyrimidines (67) (See Figure 93).

![](_page_41_Figure_8.jpeg)

bifunctional Zn-dithiocarbamate complex as catalyst (0.1 mol%)

![](_page_41_Picture_10.jpeg)

Figure 93: Three-component synthesis in the presence of a ferrocene-dithiocarbamate zinc(II) complex as the catalyst.

EtOH, 80 °C, 4 h

#### p-Toluenesulfonic Acid as the Catalytic Entity

Güngör prepared ten imidazo[1,2-a]pyrimidines with imidazole substituents in which imidazo[1,2a]pyrimidine-2-carbaldehyde, benzil, a selection of primary amines, and ammonium acetate with ptoluenesulfonic acid as catalyst were reacted in a one-pot and multicomponent reaction, under microwave irradiation. The reaction yields are medium to good, carrying many different substituents on the skeleton and the reaction conditions are optimum. The solvent utilized was ethanol, which can be considered as a green solvent. (68) (See Figures 94 and 95).

![](_page_42_Figure_2.jpeg)

**Figure 94:** Three-component reaction yields a main product and a by-product.

![](_page_42_Figure_4.jpeg)

Figure 95: A second reaction used a two-component reaction, then the third molecule gets in.

#### MISCELLANEOUS USES

#### **Proton Sensitive Organic Fluorescent Dyes**

Aydıner and Seferoglu devised a series of proton sensitive functional organic fluorescent dyes based on coumarin and imidazo[1,2-*a*]pyrimidine. They reported the synthetic effort, photophysics in different solvents, and protonation study with trifluoroacetic acid as the titrant. 7dialkylaminocoumarin-based compounds had the best fluorescent performance, and the compound having morpholine moiety was considered for its potential of colorimetric and luminescent pH sensor compound (69) (See Figures 96 and 97).

![](_page_42_Figure_10.jpeg)

![](_page_42_Figure_11.jpeg)

![](_page_43_Figure_2.jpeg)

**Figure 97:** Incorporation of an aromatic group at the 6-position with Suzuki coupling. Six compounds were reported.

Use of Selenoester Derivatives in the Treatment of Mercury(II) Chloride In situ-generated sodium selenocarboxylates

reacted with 2-(chloromethyl)imidazo[1,2a]pyrimidine in water and ethanol. Phenyl/4-tolyl selenoesters of imidazo[1,2-a]pyrimidine were obtained in water and 4-chlorophenyl/2-thienyl selenoesters of the same heterocyclic system were only obtained in ethanol. A model compound was used for studying the behavior of selenoesters with mercury(II) chloride. Reaction with HgCl<sub>2</sub> yielded a bis(imidazopyrimidinyl)selenide,

bis(imidazopyrimidinyl)diselenide and the HgCl<sub>2</sub> compound. The selenoester derivatives could be considered in the treatment of mercury(II) chloride-induced toxicity (70) (See Figure 98-100).

![](_page_43_Figure_8.jpeg)

Figure 98: Formation of chloromethylimidazo[1,2-a]pyrimidine scaffold.

![](_page_43_Figure_10.jpeg)

![](_page_43_Figure_11.jpeg)

![](_page_43_Figure_12.jpeg)

Figure 100 :Three separate experiments for the synthesis of mercury(II) selenide.

#### Heterocycles Having Long Fatty Chains

El-Sayed and coworkers reported condensed and non-condensed heterocycles having long fatty chains in the hope of obtaining surface active biological compounds. 2-Cyano-3-(dimethylamino)-N-octadecylacrylamide was the starting compound to obtain several pyrimidine, pyran, and similar heterocycles. Propylene oxide helped these compounds transform into nonionic surface active compounds. The high solubility of these compounds allowed them to be used in easy absorption. The authors report that they reduce the surface tension of liquids, and they have high wetting and emulsification property (71) (See Figure 101).

![](_page_44_Figure_5.jpeg)

Figure 101: An enaminonitrile compound reacts with 2-aminobenzimidazole, with the elimination of dimethylamine, forms the long-chained benzimidazolopyrimidine compound.

#### **CORROSION INHIBITION**

#### Corrosion Inhibition of C38 Steel in 1 M HCl

2-(m-Methoxyphenyl)imidazo[1,2-a]pyrimidine was tested along with another compound, in the corrosion inhibition study by Ghazoui et al. According to potentiodynamic polarization studies, it was found that the inhibitors employed behaved as mixed-type. The inhibitors were found to be adsorbed on the steel surface according to Langmuir isotherm. Corrosion rate decreased with the use of these two inhibitive compounds, and when the concentrations increased, the effect was Inhibitory remarkable. concentration of the compound was 1 millimolar and EI% was 95.7%(72) (See Figure 102).

![](_page_44_Figure_11.jpeg)

**Figure 102:** 2-(3-methoxyphenyl)imidazo[1,2-a]pyrimidine, used as a corrosion inhibitor.

#### Inhibitive Effect of 2-phenylimidazo[1,2a]pyrimidine-3-carbaldehyde

In hydrochloric acid solution, corrosive inhibition of the synthesized compound, 2-phenylimidazo[1,2*a*]pyrimidine-3-carbaldehyde, was tested on carbon steel. It is an expected consequence that increasing concentration of the inhibitor will lead to a decreased rate of corrosion. A mixed-type inhibition was discovered with the polarization studies. Langmuir's isotherm model was used to describe the adsorption behavior of the inhibitory entity. B3LYP/6-31+G(d,p) level of DFT was used to explain the inhibitive effects of the synthesized compound. Molecular dynamics simulation showed that the molecule could adsorb on alpha-iron(III) oxide (111) surface in hydrochloric acid through the oxygen and nitrogen atoms and pi-electrons. Experimental and theoretical results are quite harmonious with each other (73) (See Figure 103).

![](_page_44_Figure_16.jpeg)

**Figure 103:** 2-Phenylimidazo[1,2-a]pyrimidine-3-carbaldehyde.

#### Anodic Dissolution of Copper in Sodium Chloride Medium

The mechanism of anodic dissolution of copper in sodium chloride medium has not been discovered until now. 4-Amino-3-(phenyldiazenyl)benzo[4,5]imidazo[1,2-a]pyrimidin-2(1H)-one was tested for its efficiency as a corrosion inhibitor for copper, in sodium chloride solution at a concentration of 3.5% by weight. Aerated NaCl solution simulates a marine environment and adding the inhibitor or not adding it serves as a solution to uncovering the mechanistic part. According to potentiodynamic polarization

(PDP), the inhibitor is effective at a low concentration in the 3.5% NaCl solution. The inhibition efficiency was about 93% owing to adsorption and to electrochemical investigation. Electrochemical impedance spectroscopy (EIS) showed that the inhibitor acted with both diffusive and kinetically controlled processes. The adsorptive isotherm model fit with Langmuir's. The inhibitor can be used in marine environments like water coolers, desalination plants, power plants, and oil production. DFT and MD studies supported experimental part and explain more about the mode of adsorption on Cu surface (74) (See Figure 104).

![](_page_45_Figure_5.jpeg)

**Figure 104:** Preparation of the imidazopyrimidin-2(1H)-one compound and subsequent diazotization to obtain a representative corrosion inhibitor.

#### THEORY-POWERED EXPERIMENTAL STUDIES

#### Non-planarity of Imidazopyrimidines

2,4-Diphenyl-1,4-dihydrobenzo[4,5]imidazo[1,2*a*]pyrimidine and some other derivatives of it were researched both theoretically and experimentally, in the solvents, namely 1,4-dioxane and dimethylformamide. All of the compounds are nonplanar, which is indicative of great electronic and structural properties. One of the compounds are reactive at the highest owing to its lowest E(HOMO), E(LUMO) and DE. The electronic spectra of the compounds were studied in terms of the number of bands. In addition, TD-B3LYP/6-311G (d, p) in gas phase, 1,4-dioxane, and dimethylformamide showed a nice agreement with the spectra recorded (75) (See Figure 105).

![](_page_45_Figure_11.jpeg)

**Figure 105:** Structures synthesized in this study. R<sup>1</sup> is hydrogen, fluorine, bromine, methyl, nitro, or methoxy, and R<sup>2</sup> is hydrogen, fluorine, or methoxy.

#### Studying the Tautomerism of Six Imidazopyrimide Compounds by IR, NMR, and DFT

A combination of spectral analyses like FTIR and NMR and theoretical calculations were performed for six imidazopyrimidine compounds. The compounds were synthesized with the 1-aryl-4,5-dihydro-1H-imidazol-2-amine hydrobromide and diethyl phenylmalonate, and sodium methoxide was used as a base. Tautomeric transitions in the structures were investigated in the gas, solution, and crystalline states. NMR and FTIR analyses did not

help much about identification of the tautomeric forms present in the solution and also in the solid state. X-ray analysis, however, identified that all compounds were present in the keto-keto tautomeric form in the crystalline state. Theoretical studies at DFT/B3LYP/6-311++G(d,p) level was successful to indicate that, in the gas phase and in the solution, two tautomeric forms coexist. Theoretical studies and FTIR-NMR spectra allowed for the calculation of all tautomeric forms and, to a limited extent, the identification is possible (76) (See Figure 106).

![](_page_46_Figure_2.jpeg)

**Figure 106:** Synthesis of a reduced form of an imidazopyrimidine and illustration of existing tautomeric forms. Six compounds with differing R radicals were reported.

**Single-Crystal and Theoretical Study of an Imidazopyrimidine and its Thiazole Analog** Dylong and coworkers prepared imidazo[1,2*a*]pyrimidin-2-ylacetic acid and its thiazole analog, and structurally characterized with single-crystal Xray diffraction and computational studies in the form of Hirshfeld surfaces, which helped understand the intermolecular interactions of the crystals. FTIR and Raman spectra of both crystals were obtained and interpretations were made with DFT calculations and potential energy distribution (PED) analyses of computed normal vibrations. When the researchers recrystallized the compounds from deuterium oxide, they observed that the carboxylic OH groups were completely deuteriated (77) (See Figure 107).

#### 1,4-Dimethyl-2-oxopyrimido[1,2a]benzimidazole Hydrate and Some Properties of it

El-Bakri colleagues synthesized and а new 1,4-dimethyl-2-oxopyrimido[1,2compound, a]benzimidazole hydrate and used NMR and FTIR, and also used single-crystal X-ray diffraction to characterize the compound. The tricyclic core was found to be almost planar with X-ray studies. The molecules stack along the a-axis direction in headto-tail fashion through p-stacking interactions involving all three rings. DFT calculations at B3LYP/6-311 G(d,p) in gas phase helped find a very good correlation between experimental and theoretical data (78) (See Figure 108).

![](_page_46_Figure_7.jpeg)

**Figure 107:** Synthesis of imidazopyrimidine ester in low yield and basic/acidic hydrolysis of the compound to obtain a carboxylic acid.

![](_page_46_Figure_9.jpeg)

Figure 108: Formation of the benzimidazolopyrimidine and subsequent methylation of the compound.

#### REFERENCES

1. Jensen MS, Hoerrner RS, Li W, Nelson DP, Javadi GJ, Dormer PG, et al. Efficient Synthesis of a GABA A a 2,3 -Selective Allosteric Modulator via a Sequential Pd-Catalyzed Cross-Coupling Approach. J Org Chem. 2005 Jul 1;70(15):6034–9. <<u>URL></u>.

2. Oganisyan AS, Noravyan A, Dzhagatspanyan I, Nazaryan I, Akopyan A. Derivatives of condensed thienopyrimidines: Synthesis and anticonvulsant and tranquilizer activity of hydrochlorides of pyrano [4', 3': 4, 5] thieno [3, 2-e] imidazo [1, 2-a] pyrimidines. Pharm Chem J. 2007;41(11):588–90.

3. Vidal A, Ferrándiz ML, Ubeda A, Acero-Alarcón A, Sepulveda-Arques J, Alcaraz MJ. Effect of imidazo[1,2-

a]pyrimidine derivatives on leukocyte function. Inflamm Res. 2001 Jun;50(6):317–20.  $\leq$ URL>.

4. Blackaby WP, Atack JR, Bromidge F, Castro JL, Goodacre SC, Hallett DJ, et al. Imidazo[1,2-a]pyrimidines as functionally selective GABAA ligands. Bioorg Med Chem Lett. 2006 Mar;16(5):1175–9. <<u>URL></u>.

5. Jennings ASR, Lewis RT, Russell MGN, Hallett DJ, Street LJ, Castro JL, et al. Imidazo[1,2-b][1,2,4]triazines as a2/a3 subtype selective GABAA agonists for the treatment of anxiety. Bioorg Med Chem Lett. 2006 Mar;16(6):1477–80. <<u>URL></u>.

6. Zhou JP, Ding YW, Zhang HB, Xu L, Dai Y. Synthesis and anti-inflammatory activity of imidazo[1,2-a]pyrimidine derivatives. Chin Chem Lett. 2008 Jun;19(6):669–72. <<u>URL></u>.

7. Al-Tel TH, Al-Qawasmeh RA. Post Groebke– Blackburn multicomponent protocol: Synthesis of new polyfunctional imidazo[1,2-a]pyridine and imidazo[1,2a]pyrimidine derivatives as potential antimicrobial agents. Eur J Med Chem. 2010 Dec;45(12):5848–55. <u><URL></u>.

8. Fu R geng, You Q dong, Yang L, Wu W tong, Jiang C, Xu X li. Design, synthesis and bioevaluation of dihydropyrazolo[3,4-b]pyridine and benzo[4,5]imidazo[1,2-a]pyrimidine compounds as dual KSP and Aurora-A kinase inhibitors for anti-cancer agents. Bioorg Med Chem. 2010 Nov 15;18(22):8035–43. <<u>URL></u>.

9. Cosimelli B, Laneri S, Ostacolo C, Sacchi A, Severi E, Porcù E, et al. Synthesis and biological evaluation of imidazo[1,2-a]pyrimidines and imidazo[1,2-a]pyridines as new inhibitors of the Wnt/ $\beta$ -catenin signaling. Eur J Med Chem. 2014 Aug;83:45–56. <u><URL></u>.

10. Aeluri R, Alla M, Polepalli S, Jain N. Synthesis and antiproliferative activity of imidazo[1,2-a]pyrimidine Mannich bases. Eur J Med Chem. 2015 Jul;100:18–23. <<u>URL></u>.

11. Kaur N, Kaur K, Raj T, Kaur G, Singh A, Aree T, et al. One-pot synthesis of tricyclic dihydropyrimidine derivatives and their biological evaluation. Tetrahedron. 2015 Jan;71(2):332–7.  $\leq$ URL $\geq$ .

12. Velázquez-Olvera S, Salgado-Zamora H, Jiménez-Cardoso E, Campos-Aldrete ME, Pérez-González C, Ben Hadda T. In vitro anti-Giardia lamblia activity of 2-aryl-3hydroxymethyl imidazo[1,2-a]pyridines and -pyrimidines, individually and in combination with albendazole. Acta Trop. 2016 Mar;155:6-10. <<u>URL></u>.

13. Al-Lami N, Amer Z, Ali RA. Pryparation, Characterization And Biological Activity of New Derivatives of 2-Biphenyl-3-Aminomethylimidazo (1, 2-a) Pyrimidine. J Pharm Sci Res. 2018;10(12):3344–50.

14. Moszczyński-Pętkowski R, Majer J, Borkowska M, Bojarski Ł, Janowska S, Matłoka M, et al. Synthesis and characterization of novel classes of PDE10A inhibitors - 1H-1,3-benzodiazoles and imidazo[1,2-a]pyrimidines. Eur J Med Chem. 2018 Jul;155:96–116. <u><URL></u>.

15. Farag AM, Fahim AM. Synthesis, biological evaluation and DFT calculation of novel pyrazole and pyrimidine derivatives. J Mol Struct. 2019 Mar;1179:304–14. <<u>URL></u>.

16. Mantipally M, Gangireddy MR, Gundla R, Badavath VN, Mandha SR, Maddipati VC. Rational design, molecular docking and synthesis of novel homopiperazine linked imidazo[1,2-a]pyrimidine derivatives as potent cytotoxic and antimicrobial agents. Bioorg Med Chem Lett. 2019 Aug;29(16):2248–53. <<u>URL></u>.

17. Tawfik HO, El-Moselhy TF, El-Din NS, El-Hamamsy MH. Design, synthesis, and bioactivity of dihydropyrimidine derivatives as kinesin spindle protein inhibitors. Bioorg Med Chem. 2019 Dec;27(23):115126. <<u>URL></u>.

18. Elenich OV, Lytvyn RZ, Blinder OV, Skripskaya OV, Lyavinets OS, Pitkovych KhE, et al. Synthesis and Antimicrobial Activity of 3-Phenyl-1-Methylquinolin-2-One Derivatives. Pharm Chem J. 2019 Mar 15;52(12):969-74. <URL>.

19. Jismy B, Akssira M, Knez D, Guillaumet G, Gobec S, Abarbri M. Efficient synthesis and preliminary biological evaluations of trifluoromethylated imidazo[1,2- a ]pyrimidines and benzimidazo[1,2- a ]pyrimidines. New J Chem. 2019;43(25):9961–8. <<u>URL></u>.

20. D'Aria F, D'Amore VM, Di Leva FS, Amato J, Caterino M, Russomanno P, et al. Targeting the KRAS oncogene: Synthesis, physicochemical and biological evaluation of novel G-Quadruplex DNA binders. Eur J Pharm Sci. 2020 Jun;149:105337. <<u>URL></u>.

21. Zhuang J, Ma S. Recent Development of Pyrimidine-Containing Antimicrobial Agents. ChemMedChem. 2020 Oct 19;15(20):1875–86. <<u>URL></u>.

22. Zadeh EE, Vahabpour R, Beirami AD, Hajimahdi Z, Zarghi A. Novel 4-Oxo-4,10dihydrobenzo[4,5]imidazo[1,2-a]pyrimidine-3-carboxylic Acid Derivatives as HIV-1 Integrase Inhibitors: Synthesis, Docking Studies, Molecular Dynamics Simulation and Biological Activities. Med Chem. 2021 Nov;17(9):1060-71. <URL>.

23. Lima MLSO, Braga CB, Becher TB, Odriozola-Gimeno M, Torrent-Sucarrat M, Rivilla I, et al. Fluorescent Imidazo[1,2- a ]pyrimidine Compounds as Biocompatible Organic Photosensitizers that Generate Singlet Oxygen: A Potential Tool for Phototheranostics. Chem – Eur J. 2021 Apr 7;27(20):6213–22.  $\leq$  URL $\geq$ .

24. Saldabol NO, Popelis J, Slavinskaya VA. Formylation of furyl-substituted imidazo[1,2-a]pyridine, imidazo[1,2-a]pyrimidine and imidazo[2,1-b]thiazole. Chem Heterocycl Compd. 2001;37(8):1021–4. <u>URL></u>.

25. Majcen-Le Marechal A, Le Grel P, Robert A, Biškup J, Ferk V, Toplak R. Reactivity of 3-halo-2-oxopropanamides and 3-halo-2-cyano-2-hydroxypropanoates: synthesis of S and N containing heterocycles. ARKIVOC. 2001;119:132.

26. Kochergin PM, Mazur IA, Rogul'chenko GK, Aleksandrova EV, Mandrichenko BE. Investigations of Imidazoles. 101. Synthesis of 2,3-Dihydroimidazo[1,2-a]pyrimidine from 2-Aminopyrimidine. Chem Heterocycl Compd. 2001;37(2):224-6. <<u>URL></u>.

27. Chowdhury AZMS, Shibata Y. Synthesis of fused pyrimidines with N -, S -heterocyclic moieties by doubleannelation reaction. J Heterocycl Chem. 2001 May;38(3):743–7.  $\leq$  URL>. 28. Elmaati TMA. Reactions with Heterocyclic Amidines: New Routes for the Benzo [4, 5] imidazo [1, 2-a] pyrimidine, Some Pyridine, Pyran and Pyrazole Derivatives Containing the Antipyrine Moiety. Acta Chim Slov. 2002;49:721–32.

29. Katritzky AR, Xu YJ, Tu H. Regiospecific Synthesis of 3-Substituted Imidazo[1,2-a]pyridines, Imidazo[1,2-a]pyrimidines, and Imidazo[1,2-c]pyrimidine. J Org Chem. 2003 Jun 1;68(12):4935–7. <u><URL></u>.

30. El Kazzouli S, Berteina-Raboin S, Mouaddib A, Guillaumet G. Solid-phase synthesis of imidazo[1,2-a]pyridines and imidazo[1,2-a]pyrimidines. Tetrahedron Lett. 2003 Aug;44(33):6265–7. <<u>URL></u>.

31. Bari A, Feist H, Michalik D, Michalik M, Peseke K. Synthesis of Heterocyclic Compounds by Ring Transformations of 2-Formyl Pentose Glycals. Synthesis. 2004 Oct 7;2004(17):2863–8. <<u>URL></u>.

32. Sączewski J, Brzozowski Z, Gdaniec M. Synthesis of 2,3-dihydroimidazo[1,2-a]pyrimidin-5(1H)-ones by the domino Michael addition retro-ene reaction of 2-alkoxyiminoimidazolidines and acetylene carboxylates. Tetrahedron. 2005 May;61(22):5303–9. <<u>URL></u>.

33. Laquintana V, Denora N, Lopedota A, Suzuki H, Sawada M, Serra M, et al. N -Benzyl-2-(6,8-dichloro-2-(4-chlorophenyl)imidazo[1,2- a ]pyridin-3-yl)- N -(6- (7-nitrobenzo[ c ][1,2,5]oxadiazol-4-ylamino)hexyl)acetamide as a New Fluorescent Probe for Peripheral Benzodiazepine Receptor and Microglial Cell Visualization. Bioconjug Chem. 2007 Sep 1;18(5):1397–407.  $\leq$ URL>.

34. El Otmani B, El Mahi M, Bouhfid R, Essassi EM, Rohand T, Dehaen W, et al. Synthesis of new derivatives of 2, 3-dicyano-imidazo [1, 2-a] pyrimidine from 4-hydroxy-6-methylpyran-2-ones. Arkivoc. 2008;2:59–70.

35. Muzalevskiy V, Nenajdenko V, Shastin A, Balenkova E, Haufe G. Synthetic Approach to Alkoxy- $\beta$ -(trifluoromethyl)styrenes and Their Application in the Synthesis of New Trifluoromethylated Heterocycles. Synthesis. 2009 Jul;2009(13):2249–59. <<u>URL></u>.

36. Ermolat'ev DS, Svidritsky EP, Babaev EV, Van der Eycken E. Microwave-assisted synthesis of substituted 2amino-1H-imidazoles from imidazo[1,2-a]pyrimidines. Tetrahedron Lett. 2009 Sep;50(37):5218–20. <u><URL></u>.

37. Bakherad M, Keivanloo A, Kalantar Z, Jajarmi S. Pd/C-catalyzed heterocyclization during copper-free Sonogashira coupling: synthesis of 2-benzylimidazo[1,2-a]pyrimidines in water. Tetrahedron Lett. 2011 Jan;52(2):228–30. <<u>URL></u>.

38. Li Z, Hong J, Zhou X. An efficient and clean CuIcatalyzed chalcogenylation of aromatic azaheterocycles with dichalcogenides. Tetrahedron. 2011 May;67(20):3690–7. <u>< URL></u>.

39. Hu L, Zhan Z, Lei M, Hu L. Facile and Green Method for the Synthesis of 4-Amino-1,2-Dihydrobenzo [4,5]Imidazo[1,2- a ]Pyrimidine-3-Carbonitriles Catalysed by Ammonium Acetate. J Chem Res. 2012 Dec;36(12):738–9. <<u>URL></u>.

40. Borisov A, Tolmachev A, Zavada O, Kovalenko S. Synthesis of imidazo [1, 2-a] pyrazine and imidazo [1, 2-a] pyrimidine derivatives. Chem Heterocycl Compd. 2013;49(5):704–11.

41. Mehranpour AM, Zahiri M. Synthesis and characterization of new benzimidazole derivatives using 2-substituted 1,3-bis(dimethylamino)-trimethinium salts. Tetrahedron Lett. 2014 Jul;55(29):3969–71. <u><URL></u>.

42. Burbulienė MM, Dargytė J, Vainilavičius P. Interaction of methyl (2-methylsulfanyl-6-phenyl-4thioxopyrimidin-3 (4H)-yl) acetate with hydrazine hydrate. Chemija. 2014;25(1).

43. Aziz J, Piguel S. An Update on Direct C–H Bond Functionalization of Nitrogen-Containing Fused Heterocycles. Synthesis. 2017 Oct;49(20):4562–85. <URL>.

44. Bhagat SB, Telvekar VN. NBS mediated protocol for the synthesis of N -bridged fused heterocycles in water. Tetrahedron Lett. 2017 Sep;58(37):3662–6. <<u>URL></u>.

45. Rawat M, Rawat DS. Copper oxide nanoparticle catalysed synthesis of imidazo[1,2-a]pyrimidine derivatives, their optical properties and selective fluorescent sensor towards zinc ion. Tetrahedron Lett. 2018 Jun;59(24):2341–6.  $\leq$  URL>.

46. Bettanin L, Saba S, Doerner CV, Franco MS, Godoi M, Rafique J, et al. NH4I-catalyzed chalcogen(S/Se)-functionalization of 5-membered N-heteroaryls under metal-free conditions. Tetrahedron. 2018 Jul;74(29):3971–80. <<u>URL></u>.

47. Yagodinets PI, Rusnak OV, Lytvyn RZ, Skrypska OV, Pitkovych KhYe, Obushak MD. 4-(4-Acethylphenyl)-3hydroxycoumarin in the Synthesis of Nitrogen-containing Heterocycles with a Neoflavonoid Moiety. Russ J Org Chem. 2019 Aug;55(8):1145–52. <<u>URL></u>.

48. Makra Z, Puskás LG, Kanizsai I. A convenient approach for the preparation of imidazo[1,2- a ]-fused bicyclic frameworks via IBX/NIS promoted oxidative annulation. Org Biomol Chem. 2019;17(40):9001–7. <URL>.

49. Maji PK. Synthesis of Pyrimidine-Annulated Five-Membered Heterocycles: An Overview. Curr Org Chem. 2019 Dec 24;23(20):2204–69. <u><URL></u>.

50. White NA, Clagg K, Sirois LE, Lim NK, Nack WA, O'Shea PD, et al. Practical Synthesis of a Stable Precursor for Positron Emission Tomography Imaging Agent 18 F-GTP1. Org Process Res Dev. 2020 Sep 18;24(9):1690–9. <<u>URL></u>.

51. Obah Kosso AR, Kabri Y, Broggi J, Redon S, Vanelle P. Sequential Regioselective Diorganochalcogenations of Imidazo[1,2- a ]pyrimidines Using I 2 /H 3 PO 4 in Dimethylsulfoxide. J Org Chem. 2020 Mar 6;85(5):3071–81. <u><URL></u>.

52. Neto JSS, Balaguez RA, Franco MS, de Sá Machado VC, Saba S, Rafique J, et al. Trihaloisocyanuric acids in ethanol: an eco-friendly system for the regioselective halogenation of imidazo-heteroarenes. Green Chem. 2020;22(11):3410–5.  $\leq$ URL>.

53. Patel OPS, Nandwana NK, Legoabe LJ, Das BC, Kumar A. Recent Advances in Radical C–H Bond Functionalization of Imidazoheterocycles. Adv Synth Catal. 2020 Oct 21;362(20):4226–55. <u><URL></u>.

54. Fedotov VV, Ulomsky EN, Belskaya NP, Eltyshev AK, Savateev KV, Voinkov EK, et al. Benzimidazoazapurines: Design, Synthesis, and J Org Photophysical Study. Chem. 2021 Jun 18;86(12):8319-32. <u><URL></u>.

55. Kiselyov AS, Smith L. Novel one pot synthesis of polysubstituted pyrazolo[1,5-a]- and imidazo[1,2-a]pyrimidines. Tetrahedron Lett. 2006 Apr;47(15):2611–4. <<u>URL></u>.

56. Parchinsky VZ, Shuvalova O, Ushakova O, Kravchenko DV, Krasavin M. Multi-component reactions between 2-aminopyrimidine, aldehydes and isonitriles: the use of a nonpolar solvent suppresses formation of multiple products. Tetrahedron Lett. 2006 Feb;47(6):947–51. <<u>URL></u>.

57. Adib M, Sheibani E, Bijanzadeh HR, Zhu LG. A new, one-pot, multi-component synthesis of imines of 3-amino-2-arylimidazo[1,2-a]pyridines, 3-amino-2-arylimidazo[1,2-a]pyrazines, and 3-amino-2-arylimidazo[1,2-a]pyrimidines. Tetrahedron. 2008 Nov;64(47):10681–6. <<u>URL></u>.

58. Reddy MV, Byeon KR, Park SH, Kim DW. Polyethylene glycol methacrylate-grafted dicationic imidazolium-based ionic liquid: Heterogeneous catalyst for the synthesis of aryl-benzo[4,5]imidazo[1,2- a ]pyrimidine amines under solvent-free conditions. Tetrahedron. 2017 Aug;73(35):5289–96. <<u>URL></u>.

59. Hamidinasab M, Mobinikhaledi A. Organoacid-Decorated NiFe 2 O 4 Nanoparticles: An Efficient Catalyst for Green Synthesis of 2 H -Indazolo[2, 1-b]phthalazine-Triones and Pyrimido[1,2-a]benzimidazoles. ChemistrySelect. 2019 Jan 10;4(1):17–23. <<u>URL></u>.

60. Basyouni WM, Abbas SY, Abdelazeem NM, El-Bayouki KAM, El-Kady MY. Silica sulfuric acid/ethylene glycol as an efficient catalyst for the synthesis of benzo[4,5]imidazo[1,2- a ]pyrimidine-3-carbonitrile derivatives. Synth Commun. 2019 Nov 17;49(22):3112– 20. <<u>URL></u>.

61. Damghani FK, Pourmousavi SA, Kiyani H. Sulfonic Acid-Functionalized Magnetic Nanoparticles as an Efficient Catalyst for the Synthesis of Benzo[4, 5]imidazo[1, 2a]pyrimidine Derivatives, 2-Aminobenzothia Zolomethylnaphthols and 1-Amidoalkyl-2-naphthols. Curr Org Synth. 2019 Dec 26;16(7):1040–54. <u><URL></u>.

62. Verma P, Pal S, Chauhan S, Mishra A, Sinha I, Singh S, et al. Starch functionalized magnetite nanoparticles: A green, biocatalyst for one-pot multicomponent synthesis of imidazopyrimidine derivatives in aqueous medium under ultrasound irradiation. J Mol Struct. 2020 Mar;1203:127410. <<u>URL></u>.

63. Güngör T. One pot, multicomponent protocol for the synthesis of novel imidazo[1,2-a]pyrimidine-based pyran analogs: a potential biological scaffold. Monatshefte Für Chem - Chem Mon. 2020 May;151(5):781–9. <<u>URL></u>.

64. Akbari A, Dekamin MG, Yaghoubi A, Naimi-Jamal propylsulfonic MR Novel magnetic acid-anchored isocyanurate-based periodic mesoporous organosilica (Iron oxide@PMO-ICS-PrSO3H) as a highly efficient and reusable synthesis nanoreactor for the sustainable of imidazopyrimidine derivatives. 2020 Sci Rep. Dec;10(1):10646. </

65. Chavan P, Bangale S, Pansare D, Shelke R, Jadhav S, Tupare S, et al. Synthesis of substituted pyrimidine using ZnFe 2 O 4 nanocatalyst via one pot multi-component reaction ultrasonic irradiation. J Heterocycl Chem. 2020 Jul 11; jhet. 4048. <<u>URL></u>.

66. Bakhshali-Dehkordi R, Ghasemzadeh MA. Fe3O4@TiO2@ILs-ZIF-8 Nanocomposite: A Robust Catalyst for the Synthesis of Benzo[4,5]imidazo[1,2a]pyrimidines. J Mol Struct. 2021 Jul;1236:130298. <URL>.

67. Anamika, Yadav CL, Drew MGB, Kumar K, Singh N. Ferrocene-Functionalized Dithiocarbamate Zinc(II) Complexes as Efficient Bifunctional Catalysts for the One-Pot Synthesis of Chromene and Imidazopyrimidine Derivatives via Knoevenagel Condensation Reaction. Inorg Chem. 2021 May 3;60(9):6446–62. <<u>URL></u>.

68. Güngör T. Microwave assisted, sequential twostep, one-pot synthesis of novel imidazo [1, 2-a] pyrimidine containing tri/tetrasubstituted imidazole derivatives. Turk J Chem. 2021;45(1):219–30. <u><URL></u>.

69. Aydıner B, Seferoğlu Z. Proton Sensitive Functional Organic Fluorescent Dyes Based on Coumarinimidazo[1,2- a ]pyrimidine; Syntheses, Photophysical Properties, and Investigation of Protonation Ability: Proton Sensitive Functional Organic Fluorescent Dyes Based on Coumarin-imidazo[1,2-a]pyrimidine; Syntheses, Photophysical Properties, and Investigatio. Eur J Org Chem. 2018 Nov 25;2018(43):5921–34. <URL>.

70. Sharma N, Kumar S, Kumar S, Mehta SK, Bhasin KK. Synthesis and characterization of fused imidazole heterocyclic selenoesters and their application for chemical detoxification of HgCl 2. New J Chem. 2018;42(4):2702–10.  $\leq$  URL $\geq$ .

71. El-Sayed R, Katouah HA. Synthesis of Pyrimidine and Pyran Derivatives with the Related Systems and the Study of Their Behavior in the Liquid Solutions. J Heterocycl Chem. 2019 Aug;56(8):2134–43. <u><URL></u>.

72. Ghazoui A, Saddik R, Benchat N, Guenbour M, Hammouti B, Al-Deyab S, et al. Comparative study of pyridine and pyrimidine derivatives as corrosion inhibitors of C38 steel in molar HCl. Int J Electrochem Sci. 2012;7:7080–97. <u><URL></u>.

73. Ech-chihbi E, Belghiti ME, Salim R, Oudda H, Taleb M, Benchat N, et al. Experimental and computational studies on the inhibition performance of the organic compound "2-phenylimidazo [1,2-a]pyrimidine-3-carbaldehyde" against the corrosion of carbon steel in 1.0 M HCl solution. Surf Interfaces. 2017 Dec;9:206–17. <URL>.

74. Pareek S, Jain D, Hussain S, Biswas A, Shrivastava R, Parida SK, et al. A new insight into corrosion inhibition mechanism of copper in aerated 3.5 wt.% NaCl solution by eco-friendly Imidazopyrimidine Dye:

experimental and theoretical approach. Chem Eng J. 2019 Feb;358:725–42.  $\leq$ URL $\geq$ .

75. Elshakre ME, Moustafa H, Hassaneen HuwaidaME, Moussa AbdelrahimZ. Electronic spectra and DFT calculations of some pyrimido[1,2-a]benzimidazole derivatives. Spectrochim Acta A Mol Biomol Spectrosc. 2015 Jun;145:1–14. <<u>URL></u>.

76. Wysocki W, Karczmarzyk Z, Rządkowska M, Szacoń E, Matosiuk D, Urbańczyk-Lipkowska Z, et al. Structural investigations of a series of 1,6-aryl-7-hydroxy-2,3-dihydroimidazo[1,2-a]pyrimidin-5(1H)-ones with potential antinociceptive activity. J Mol Struct. 2015 Aug;1094:137–47. <<u>URL></u>.

77. Dylong A, Goldeman W, Sowa M, Ślepokura K, Drożdżewski P, Matczak-Jon E. Synthesis, crystal structures and spectral characterization of imidazo[1,2-a]pyrimidin-2-yl-acetic acid and related analog with imidazo[2,1-b]thiazole ring. J Mol Struct. 2016 Aug;1117:153–63. <<u>URL></u>.

78. El Bakri Y, Anouar EH, Ramli Y, Essassi EM, Mague JT. Synthesis, crystal structure, spectroscopic characterization, Hirshfeld surface analysis, and DFT calculations of 1,4-dimethyl-2-oxo-pyrimido[1,2-a]benzimidazole hydrate. J Mol Struct. 2018 Jan;1152:154–62. <<u>URL></u>.